

PROTECT YOUR PATIENTS



PATIENT SAFE[®] SYRINGE

The unique luer guard of the Patient Safe[®] syringe promotes aseptic technique

- Reduces the risk of luer tip contact contamination and medication contamination
- Reduces the risk of bloodstream infections resulting from catheter hub contamination
- Compatible with most available luer-fitting devices



The luer tip of the standard syringe is exposed, increasing the risk of luer tip contamination and medication contamination

EASY TO USE – NO TRAINING NECESSARY

Use Patient Safe[®] syringes the same as you would standard syringes. Available in a range of sizes. Please see our website for additional information.



RETRACTABLE
TECHNOLOGIES, INC.
www.retractable.com

M031 REV.2 EN, JUL 2016
© 2010, 2016 Retractable Technologies, Inc. All rights reserved

511 Lobo Lane · Little Elm, Texas 75068-5295 · USA
Tel: (972) 294-1770 · Fax: (972) 294-4400
Toll Free: 1-888-703-1010